

**THUASNE**

## LombaStab

fr	Ceinture de soutien lombaire – "Serrage facile et précis"	5
en	Lumbar support belt - "Easy, precise tightening"	6
de	LWS-Stabilisierungsorthese- „einfaches und präzises Justieren"	7
nl	Lumbale rugbrace - "Eenvoudige en nauwkeurige aanspanning"	8
it	Cintura di sostegno lombare – "Serraggio facile e preciso"	9
es	Faja de sujeción lumbar – "Fijación fácil y precisa"	10
pt	Cinta de suporte lombar – "Aperto fácil e preciso"	11
da	Rygstøttebælte – "Nem og præcis stramning"	12
fi	Lantiotukivö – Helposti kiristettävä"	14
sv	Stödbälte för ryggen – "enkel och precis åtdragning"	15
el	Ζώνη οσφυϊκής υποστήριξης - «Εύκολη και ακριβής σύσφιξη»	16
cs	Bederní pás – „Jednoduché a přesné stažení"	17
pl	Orteza lędźwiowo-krzyżowa – „Łatwe i precyzyjne mocowanie"	18
lv	Muguras lejasdaļas josta – „Viegla un precīza savilkšana"	19
lt	Palaikomas juosmens diržas – „Lengvai ir tiksliai suveržiamas"	20
et	Nimpepiirkonna tugivöö – „Lihtne ja täpne pingutus"	21
sl	Oporni križni pas – „Enostavno in natančno vpetje"	22
sk	Podporný driekový pás – „Jednoduché a rovnomerné utiahnutie"	23
hu	Ágyéki fűző – „Egyszerű és pontos rögzítés"	24
bg	Лумбален колан – „лесно и точно затягане"	25
ro	Centură de susținere lombară – „Strângere ușoară și precisă"	26
ru	Опорный поясничный корсет – «Легко и просто затянуть»	28
hr	Lumbalni potporni pojas – „Jednostavno i precizno stezanje"	29
zh	易于调节的加强型腰部稳定护带	30
ar	حزام لدعم أسفل الظهر - "شد سهل ودقيق"	31

	1	60 - 80 cm
	2	75 - 95 cm
	3	90 - 110 cm
	4	105 - 125 cm
	5	120 - 140 cm

<b>fr</b>	Tour de taille	Hauteurs : 21 cm - 26 cm
<b>en</b>	Waist circumference	Heights: 21 cm - 26 cm
<b>de</b>	Tailenumfang	Höhen: 21 cm - 26 cm
<b>nl</b>	Tailleomvang	Hoogtes: 21 cm - 26 cm
<b>it</b>	Circonferenza vita	Altezze: 21 cm - 26 cm
<b>es</b>	Contorno de cintura	Alturas: 21 cm - 26 cm
<b>pt</b>	Perímetro de cintura	Alturas: 21cm - 26cm
<b>da</b>	Taljemål	Højde: 21cm - 26cm
<b>fi</b>	Vyötärönympäry	Korkeudet: 21 cm - 26 cm
<b>sv</b>	Midjemått	Höjd: 21cm - 26cm
<b>el</b>	Περίμετρος μέσης	Ύψη: 21 εκ. - 26 εκ.
<b>cs</b>	Obvod pasu	Výšky: 21 cm - 26 cm
<b>pl</b>	Obwód w talii	Wysokości: 21 cm - 26 cm
<b>lv</b>	Vidukļa apkārtmērs	Augstumi: 21 cm - 26 cm
<b>lt</b>	Juosmens apimtis	Aukštis: 21 cm - 26 cm
<b>et</b>	Vööübermõõt	Kõrgused: 21 cm - 26 cm
<b>sl</b>	Obseg pasu	Višina: 21 cm - 26 cm
<b>sk</b>	Obvod pása	Výšky: 21 cm - 26 cm
<b>hu</b>	Derékbőség	Magasság: 21 cm - 26 cm
<b>bg</b>	Обиколка на талията	Височини: 21 cm - 26 cm
<b>ro</b>	Circonfereința taliei	Înălțimi: 21 cm - 26 cm
<b>ru</b>	Обхват талии	Высота: 21 см - 26 см
<b>hr</b>	Opseg struka	Visine: 21 cm - 26 cm
<b>zh</b>	腰围	高度：21 cm - 26 cm
<b>ar</b>	محيط الخصر	الارتفاعات: 21 سم - 26 سم

Compression	Maintien	Éducation posturale
Compression	Support	Postural education
Kompression	Halt	Erlernen der Körperhaltung
Compressie	Fixatie	Voorlichting over houding
Compressione	Sostegno	Educazione posturale
Compresión	Sujeción	Educación postural
Compressão	Manutenção	Educação postural
Kompression	Støtte	Postural genoptræning
Kompressio	Tuki	Asentokuntoutus
Kompression	Stöd	Postural träning
Συμπίεση	Συγκράτηση	Σωστή στάση σώματος
Komprese	Stabilizace	Správné držení těla
Kompresja	Podtrzymywanie	Korekcja postawy
Kompresija	Atbalsts	Stājas korekcija
Kompresija	Palaikymas	Laikysenos taisymas
Kompressioon	Hooldus	Rühtreening
Kompresija	Oprijem	Učenje telesne drže
Kompresia	Podpora	Cvičenie posturálnych svalov
Szorítás	Támasz	Tartásjavító oktatás
Компресия	Поддръжка	Приучаване към правилна стойка
Compressie	Sușținere	Educație posturală
Компресия	Стабилизация	Исправление осанки
Kompresija	Potpora	Učenje o držanju
压紧	穩定	姿勢联系
الضغط	الدعم	تدريب وضعي

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## en

# LUMBAR SUPPORT BELT - "EASY, PRECISE TIGHTENING"\*

## Description/purpose

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

## Composition

Elastic fabric. ①

Back panel with Coolmax lining. ②

Abdominal panels with finger-loops.

Handles. ③

Laces.

Textile components: polyamide - polyester - elastane-viscose - polyurethane - polyéthylène.

Non-textile components: brushed tempered steel - polypropylène - SEBS - polyamide - polyuréthane - polyoxyméthylène.

## Properties/Mode of action

The device relieves low back pain, supporting the lumbar spine: compartmental lumbar unloading is therefore achieved via increased intra-abdominal pressure.

This device uses back pads, implementing a pulley-like system to distribute pressure.

## Indications

- Non-specific low back pain (acute, subacute and chronic phases).
- Spinal stenosis.

## Contraindications

Do not use on pregnant women.

Do not use in the event of known allergy to any of the components.

Do not use in the event of a hiatus hernia.

Do not use in the event of bone cancer with spinal metastases.

Do not use in the event of circulatory, pulmonary or cardiovascular problems in patients for whom an increase in arterial blood pressure is not recommended.

## Precautions

Verify the product's integrity before every use.

Wear the device above thin clothing.

Do not use the product if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

Strictly comply with your healthcare professional's prescription and recommendations for use.

It is recommended that a healthcare professional adjust the stays to the shape of the patient's back before the belt is worn for the first time.

For hygiene, security and performance reasons, do not re-use the product for another patient.

Do not apply the product in direct contact with skin (risk of pinching).

For safety reasons, do not wear the belt close to machines (risk of becoming caught).

Do not wear the product in a medical imaging machine.

Do not wear the product while sleeping.

Consult a health professional if discomfort or pain occurs.

The use of a medical device by a child, should be done under the supervision of an adult or a healthcare professional.

## Undesirable side-effects

This device can cause skin reactions (redness, itching, irritations...).

Any serious incidents related to the product should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

## Instructions for use/Application

Open the belt and lay it out completely flat so that you can see the handles.

Check that the two tightening handles are fully retracted

If necessary, unhook them and separate the two back parts by pulling the ends of the belt, then fasten the handles again.

Position the belt around the waist, using the finger loops located on the front panels. ④

The tightening back pads must be positioned outwards (i.e. : not against the skin).

The lacing section must be centered over the spine: the pads should be positioned on either side of the spine. The bottom of the belt must stop at the top of the crease between the buttocks.

The Thuasne logo indicates which way up the belt should be put on.

Fasten the belt: right side over left side, without pulling too tight or applying any particular pressure. ⑤ ⑥

Adjust the tightness using the two handles: pull them forwards then fold them flat against the front of the belt, where they attach like a touch fastener. ⑦ ⑧

These two handles can be used to adjust the level of compression as required throughout the day

## Care

Close the self-fastening tabs before washing. Machine washable at 30°C (delicate programme). Do not machine wash the device more than 10 times. Do not dry clean. Do not use detergents, fabric softeners or aggressive products (products containing chlorine). Do not tumble-dry. Do not iron. Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). Dry flat.

## Storage

Store at room temperature, preferably in the original packaging.

## Disposal

Dispose of the product with household waste (no special disposal conditions).

\*Study conducted in-house on a panel of 13 people, July 2019

This medical device is a regulated health product with CE marking.

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©Thuasne - 2122102 (2022-06)
THUASNE UK Ltd**
Thuasne SAS - SIREN/RCS Nanterre 542 091 186